

RAB21 Antibody, Biotin conjugated



PACO59560

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Regulates integrin internalization and recycling, but does not influence the traffic of endosomally translocated receptors in general. As a result, may regulate cell adhesion and migration. During the mitosis of adherent cells, controls the endosomal trafficking of integrins which is required for the successful completion of cytokinesis. Involved in neurite growth.

Gene ID:

RAB21

Uniprot

Q9UL25

Synonyms:

Ras-related protein Rab-21, RAB21, KIAA0118

Immunogen:

Recombinant Human Ras-related protein Rab-21 protein (141-222AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A