

# CGN Antibody, HRP conjugated



PACO59693

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Probably plays a role in the formation and regulation of the tight junction (TJ) paracellular permeability barrier.

**Gene ID:**

CGN

**Uniprot**

Q9P2M7

**Synonyms:**

Cingulin, CGN, KIAA1319

**Immunogen:**

Recombinant Human Cingulin protein (92-252AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A