

CGN Antibody, Biotin conjugated



PACO59695

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probably plays a role in the formation and regulation of the tight junction (TJ) paracellular permeability barrier.

Gene ID:

CGN

Uniprot

Q9P2M7

Synonyms:

Cingulin, CGN, KIAA1319

Immunogen:

Recombinant Human Cingulin protein (92-252AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A