

DHCR7 Antibody, Biotin conjugated



PACO59719

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Production of cholesterol by reduction of C7-C8 double bond of 7-dehydrocholesterol (7-DHC).

Gene ID:

DHCR7

Uniprot

Q9UBM7

Synonyms:

7-dehydrocholesterol reductase (7-DHC reductase) (EC 1.3.1.21) (Putative sterol reductase SR-2) (Sterol Delta(7)-reductase), DHCR7, D7SR

Immunogen:

Recombinant Human 7-dehydrocholesterol reductase protein (63-147AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A