

POU1F1 Antibody, HRP conjugated



PACO59853

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcription factor involved in the specification of the lactotrope, somatotrope, and thyrotrope phenotypes in the developing anterior pituitary. Specifically binds to the consensus sequence 5'-TAAAT-3'. Activates growth hormone and prolactin genes.

Gene ID:

POU1F1

Uniprot

P28069

Synonyms:

Pituitary-specific positive transcription factor 1 (PIT-1) (Growth hormone factor 1) (GHF-1), POU1F1, GHF1 PIT1

Immunogen:

Recombinant Human Pituitary-specific positive transcription factor 1 protein (34-130AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A