

# POU1F1 Antibody, FITC conjugated



PACO59854

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Transcription factor involved in the specification of the lactotrope, somatotrope, and thyrotrope phenotypes in the developing anterior pituitary. Specifically binds to the consensus sequence 5'-TAAAT-3'. Activates growth hormone and prolactin genes.

**Gene ID:**

POU1F1

**Uniprot**

P28069

**Synonyms:**

Pituitary-specific positive transcription factor 1 (PIT-1) (Growth hormone factor 1) (GHF-1), POU1F1, GHF1 PIT1

**Immunogen:**

Recombinant Human Pituitary-specific positive transcription factor 1 protein (34-130AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A