

# SCLY Antibody, FITC conjugated



PACO59882

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Catalyzes the decomposition of L-selenocysteine to L-alanine and elemental selenium.

**Gene ID:**

SCLY

**Uniprot**

Q96I15

**Synonyms:**

Selenocysteine lyase (hSCL) (EC 4.4.1.16), SCLY, SCL

**Immunogen:**

Recombinant Human Selenocysteine lyase protein (3-138AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A