

# TBX4 Antibody



PACO59900

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

Involved in the transcriptional regulation of genes required for mesoderm differentiation. Probably plays a role in limb pattern formation.

**Gene ID:**

TBX4

**Uniprot**

P57082

**Synonyms:**

T-box transcription factor TBX4 (T-box protein 4), TBX4

**Immunogen:**


Recombinant Human T-box transcription factor TBX4 protein (355-511AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in: A549 whole cell lysate, HepG2 whole cell lysate. All lanes: TBX4 antibody at 4.3 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 61 kDa. Observed band size: 61 kDa.