

UMODL1 Antibody, FITC conjugated



PACO59934

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	UMODL1
Source:	Uniprot
Rabbit	Q5DID0
Isotype:	Synonyms:
IgG	Uromodulin-like 1 (Olfactorin), UMODL1
Applications:	Immunogen:
ELISA	Recombinant Human Uromodulin-like 1 protein (537-654AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A