

UMODL1 Antibody, Biotin conjugated



PACO59935

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

UMODL1

Uniprot

Q5DID0

Synonyms:

Uromodulin-like 1 (Olfactorin), UMODL1

Immunogen:

Recombinant Human Uromodulin-like 1 protein (537-654AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A