

PTPRK Antibody, Biotin conjugated



PACO59987

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Regulation of processes involving cell contact and adhesion such as growth control, tumor invasion, and metastasis. Negative regulator of EGFR signaling pathway. Forms complexes with beta-catenin and gamma-catenin/plakoglobin. Beta-catenin may be a substrate for the catalytic activity of PTPRK/PTP-kappa.

Gene ID:

PTPRK

Uniprot

Q15262

Synonyms:

Receptor-type tyrosine-protein phosphatase κ (Protein-tyrosine phosphatase kappa) (R-PTP-kappa) (EC 3.1.3.48), PTPRK, PTPK

Immunogen:

Recombinant Human Receptor-type tyrosine-protein phosphatase κ protein (589-749AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A