

# IQGAP2 Antibody, Biotin conjugated



PACO60003

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Binds to activated CDC42 and RAC1 but does not seem to stimulate their GTPase activity. Associates with calmodulin.

**Gene ID:**

IQGAP2

**Uniprot**

Q13576

**Synonyms:**

Ras GTPase-activating-like protein IQGAP2, IQGAP2

**Immunogen:**

Recombinant Human Ras GTPase-activating-like protein IQGAP2 protein (1245-1395AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A