

# DAZAP2 Antibody, Biotin conjugated



PACO60031

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

nuclear speck, identical protein binding, WW domain binding.

**Gene ID:**

DAZAP2

**Uniprot**

Q15038

**Synonyms:**

DAZ-associated protein 2 (Deleted in azoospermia-associated protein 2), DAZAP2, KIAA0058

**Immunogen:**

Recombinant Human DAZ-associated protein 2 protein (1-168AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A