

INF2 Antibody, HRP conjugated



PACO60049

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Severs actin filaments and accelerates their polymerization and depolymerization.

Gene ID:

INF2

Uniprot

Q27J81

Synonyms:

Inverted formin-2 (HBEBP2-binding protein C), INF2, C14orf151 C14orf173

Immunogen:

Recombinant Human Inverted formin-2 protein (926-1047AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A