

# HYI Antibody, Biotin conjugated



PACO60087

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Catalyzes the reversible isomerization between hydroxypyruvate and 2-hydroxy-3-oxopropanoate (also termed tartronate semialdehyde).

**Gene ID:**

HYI

**Uniprot**

Q5T013

**Synonyms:**

Putative hydroxypyruvate isomerase (EC 5.3.1.22) (Endothelial cell apoptosis protein E-CE1), HYI

**Immunogen:**

Recombinant Human Putative hydroxypyruvate isomerase protein (1-174AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A