

VARS2 Antibody, FITC conjugated



PACO60106

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

valine-tRNA ligase activity, valyl-tRNA aminoacylation.

Gene ID:

VARS2

Uniprot

Q5ST30

Synonyms:

Valine--tRNA ligase, mitochondrial (EC 6.1.1.9) (Valyl-tRNA synthetase) (ValRS) (Valyl-tRNA synthetase-like), VARS2, KIAA1885 VARS2L VARSL

Immunogen:

Recombinant Human Valine--tRNA ligase, mitochondrial protein (830-981AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A