

# THEMIS2 Antibody, HRP conjugated



PACO60113

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May constitute a control point in macrophage inflammatory response, promoting LPS-induced TNF production.

**Gene ID:**

THEMIS2

**Uniprot**

Q5TEJ8

**Synonyms:**

Protein THEMIS2 (Induced by contact to basement membrane 1 protein) (Protein ICB-1) (Thymocyte-expressed molecule involved in selection protein 2), THEMIS2, C1orf38 ICB1

**Immunogen:**

Recombinant Human Protein THEMIS2 protein (62-177AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A