

THEMIS2 Antibody, Biotin conjugated



PACO60115

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May constitute a control point in macrophage inflammatory response, promoting LPS-induced TNF production.

Gene ID:

THEMIS2

Uniprot

Q5TEJ8

Synonyms:

Protein THEMIS2 (Induced by contact to basement membrane 1 protein) (Protein ICB-1) (Thymocyte-expressed molecule involved in selection protein 2), THEMIS2, C1orf38 ICB1

Immunogen:

Recombinant Human Protein THEMIS2 protein (62-177AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A