

POF1B Antibody



PACO60256

Product Information

Size:

50ug

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Plays a key role in the organization of epithelial monolayers by regulating the actin cytoskeleton. May be involved in ovary development.

Gene ID:

POF1B

Uniprot

Q8WVV4

Synonyms:

Protein POF1B (Premature ovarian failure protein 1B), POF1B


Immunogen:

Recombinant Human Protein POF1B protein (532-589AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: K562 whole cell lysate, Rat lung tissue. All lanes: POF1B antibody at 3.4 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 69, 35 kDa. Observed band size: 69 kDa.