

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Regulates, cooperatively with NONO and SFPQ, androgen receptor-mediated gene transcription activity in Sertoli cell line. Binds to poly(A), poly(G) and poly(U) RNA homopolymers. Regulates the circadian clock by repressing the transcriptional activator activity of the CLOCK-ARNTL/BMAL1 heterodimer. Together with NONO, required for the formation of nuclear paraspeckles.

**Gene ID:**

PSPC1

**Uniprot**

Q8WXF1

**Synonyms:**

Paraspeckle component 1 (Paraspeckle protein 1), PSPC1, PSP1

**Immunogen:**

Recombinant Human Paraspeckle component 1 protein (191-299AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A