

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

IRBP shuttles 11-cis and all trans retinoids between the retinol isomerase in the pigment epithelium and the visual pigments in the photoreceptor cells of the retina.

Gene ID:

RBP3

Uniprot

P10745

Synonyms:

Retinol-binding protein 3 (Interphotoreceptor retinoid-binding protein) (IRBP) (Interstitial retinol-binding protein), RBP3

Immunogen:

Recombinant Human Retinol-binding protein 3 protein (1146-1246AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A