

rbm47 Antibody



PACO60913

Product Information

Size:

50ug

Reactivity:

Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

head development, regulation of anterior head development.

Gene ID:

rbm47

Uniprot

A2BIC5

Synonyms:

RNA-binding protein 47 (RNA-binding motif protein 47), Rbm47

Immunogen:

Recombinant Zebrafish RNA-binding motif protein 47 protein (390-524AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A