

PACO60917

---

## Product Information

**Size:**

50ug

**Reactivity:**

Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

May be involved in early neuronal development.

**Gene ID:**

mturn

**Uniprot**

P0C8M4

**Synonyms:**

Maturin (Maturin neural progenitor differentiation regulator protein homolog), mturn

**Immunogen:**


Recombinant Zebrafish Maturin protein (1-133AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in Recombinant protein. All lanes: mturn antibody at 2 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 46 kDa. Observed band size: 50 kDa.