

NSUN5 Antibody, FITC conjugated



PACO61428

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

S-adenosyl-L-methionine-dependent methyltransferase that specifically methylates the C(5) position of cytosine 3782 in 28S rRNA.

Gene ID:

NSUN5

Uniprot

Q96P11

Synonyms:

Probable 28S rRNA (cytosine-C(5))-methyltransferase, NOL1-related protein, NOL1R, NOL1/NOP2/Sun domain family member 5, Williams-Beuren syndrome chromosomal region 20A protein, NSUN5, NSUN5A, WBSCR20, WBSCR20A

Immunogen:

Recombinant Human Probable 28S rRNA (cytosine-C(5))-methyltransferase protein (362-466AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A