

TSPEAR Antibody



PACO61434

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

May play a role in development or function of the auditory system.

Gene ID:

TSPEAR

Uniprot

Q8WU66

Synonyms:

Thrombospondin-type laminin G domain and EAR repeat-containing protein, TSP-EAR, TSPEAR, C21orf29


Immunogen:

Recombinant Human Thrombospondin-type laminin G domain and EAR repeat-containing protein (248-433AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: Raji whole cell lysate. All lanes: TSPEAR antibody at 3.5 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 75, 67 kDa. Observed band size: 75 kDa.