

YY1AP1 Antibody, FITC conjugated



PACO61464

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Associates with the INO80 chromatin remodeling complex, which is responsible for transcriptional regulation, DNA repair, and replication. Enhances transcription activation by YY1. Plays a role in cell cycle regulation.

Gene ID:

YY1AP1

Uniprot

Q9H869

Synonyms:

YY1-associated protein 1, Hepatocellular carcinoma susceptibility protein1, YY1AP1

Immunogen:

Recombinant Human YY1-associated protein 1 protein (335-466AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A