

NRL Antibody, Biotin conjugated



PACO61577

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcription factor which regulates the expression of several rod-specific genes, including RHO and PDE6B.

Gene ID:

NRL

Uniprot

P54845

Synonyms:

Neural retina-specific leucine zipper protein, NRL, NRL, D14S46E

Immunogen:

Recombinant Human Neural retina-specific leucine zipper protein (24-124AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A