

# SPINK5 Antibody, FITC conjugated



PACO61688

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Serine protease inhibitor, probably important for the anti-inflammatory and/or antimicrobial protection of mucous epithelia. Contribute to the integrity and protective barrier function of the skin by regulating the activity of defense-activating and desquamation-involved proteases. Inhibits KLK5, it's major target, in a pH-dependent manner. Inhibits KLK7, KLK14 CASP14, and trypsin.

**Gene ID:**

SPINK5

**Uniprot**

Q9NQ38

**Synonyms:**

Serine protease inhibitor Kazal-type 5, Lympho-epithelial Kazal-type-related inhibitor, LEKTI, , Hemofiltrate peptide HF6478, Hemofiltrate peptide HF7665, SPINK5

**Immunogen:**

Recombinant Human Serine protease inhibitor Kazal-type 5 protein (722-847AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A