

TTC39B Antibody




PACO61714

Product Information

Size:	Protein Background:
50ul	
Reactivity:	Gene ID:
Human	TTC39B
Source:	Uniprot
Rabbit	Q5VTQ0
Isotype:	Synonyms:
IgG	Tetratricopeptide repeat protein 39B (TPR repeat protein 39B), TTC39B, C9orf52
Applications:	Immunogen:
ELISA, WB	Recombinant Human Tetratricopeptide repeat protein 39Bprotein (10-122AA).
Recommended dilutions:	Storage:
ELISA:1:2000-1:10000, WB:1:500-1:5000	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: HeLa whole cell lysate, MCF-7 whole cell lysate. All lanes: TTC39B antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 77, 67, 76, 60, 70, 23 kDa. Observed band size: 77 kDa.