

# ORF7 Antibody



PACO61746

---

## Product Information

**Size:**

50ul

**Reactivity:**

Porcine reproductive and respiratory syndrome virus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

**Gene ID:**

ORF7

**Uniprot**

V5N6H2

**Synonyms:**

ORF7

**Immunogen:**


Recombinant Porcine reproductive and respiratory syndrome virus Nucleocapsid protein (1-123AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in Recombinant protein. All lanes: ORF7 antibody at 2.1 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 18 kDa. Observed band size: 18 kDa.