

PACO61766

---

## Product Information

**Size:**

50ul

**Reactivity:**

Zebrafish

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

**Gene ID:**

src

**Uniprot**

Q1JPZ3

**Synonyms:**

Proto-oncogene tyrosine-protein kinase Src (EC 2.7.10.2) (Proto-oncogene c-Src) (pp60c-src) (p60-Src), src

**Immunogen:**


Recombinant Zebrafish Proto-oncogene tyrosine-protein kinase Src protein (2-146AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in: Zebrafish tissue. All lanes: src antibody at 2 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 61 kDa. Observed band size: 61 kDa.