

PACO61774

Product Information

Size:

50ul

Reactivity:

Mycoplasma pneumoniae

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Gene ID:

pepP

Uniprot

P75313

Synonyms:

Putative Xaa-Pro aminopeptidase (X-Pro aminopeptidase) (EC 3.4.11.9)
(Aminoacylproline aminopeptidase) (Aminopeptidase P) (APP), pepP


Immunogen:

Recombinant Mycoplasma pneumoniae Putative Xaa-Pro aminopeptidase protein (1-354AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: pepP antibody at 3.8 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 58 kDa. Observed band size: 58 kDa.