

MPN_083 Antibody



PACO61778

Product Information

Size:

50ul

Reactivity:

Mycoplasma pneumoniae

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:**Gene ID:**

MPN_083

Uniprot

P75610

Synonyms:

Uncharacterized lipoprotein MPN_083


Immunogen:

Recombinant Mycoplasma pneumoniae Uncharacterized lipoprotein MPN_083 protein (22-242AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: MPN_083 at 2.5 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 43 kDa. Observed band size: 43 kDa.