

PACO61806

---

## Product Information

**Size:**

50ul

**Reactivity:**

Vibrio vulnificus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

Bacterial hemolysins are exotoxins that attack blood cell membranes and cause cell rupture by mechanisms not clearly defined.

**Gene ID:**

whA

**Uniprot**

P19247

**Synonyms:**

Cytolysin, vvhA

**Immunogen:**


Recombinant Vibrio vulnificus Cytolysin protein (21-471AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in Recombinant protein. All lanes: vvhA antibody at 2.3 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 58 kDa. Observed band size: 58 kDa.