

C Antibody



PACO61810

Product Information

Size:

50ul

Reactivity:

Rabies virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Gene ID:

C

Uniprot

P15201

Synonyms:

Putative uncharacterized protein C, C


Immunogen:

Recombinant Rabies virus Putative uncharacterized protein C protein (1-69AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: C antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 44 kDa. Observed band size: 44 kDa.