

C Antibody, Biotin conjugated



PACO61813

Product Information

Size:

50ul

Reactivity:

Rabies virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

C

Uniprot

P15201

Synonyms:

Putative uncharacterized protein C, C

Immunogen:

Recombinant Rabies virus Putative uncharacterized protein C protein (1-69AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A