

PACO61822

Product Information

Size:

50ul

Reactivity:

Vibrio cholerae serotype O1

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Increases the permeability of the small intestine mucosa by affecting the structure of intercellular tight junctions (zonula occludens).

Gene ID:

zot

Uniprot

P38442

Synonyms:

Zona occludens toxin (Zonular occludens toxin) (Zot), zot


Immunogen:

Recombinant Vibrio cholerae serotype O1 Zona occludens toxin protein (1-399AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: zot antibody at 3.5 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 52 kDa. Observed band size: 52 kDa.