

VP3 Antibody



PACO61834

Product Information

Size:

50ul

Reactivity:

Chicken anemia virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

May act as transcriptional regulator. Induces apoptosis in infected cells. Element of infectious replication cycle.

Gene ID:

VP3

Uniprot

P54095

Synonyms:

Apoptin, VP3


Immunogen:

Recombinant Chicken anemia virus Apoptin protein (1-120AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: VP3 antibody at 3.2 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 23 kDa. Observed band size: 23 kDa.