

pmpD Antibody



PACO61850

Product Information

Size:

50ul

Reactivity:

Chlamydia trachomatis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Gene ID:

pmpD

Uniprot

O84818

Synonyms:

Probable outer membrane protein PmpD (Polymorphic membrane protein D), pmpD


Immunogen:

Recombinant Chlamydia trachomatis Probable outer membrane protein PmpD protein (1244-1531AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in Recombinant protein. All lanes: pmpD antibody at 3.9 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 47 kDa. Observed band size: 47 kDa.