

PTBP3 Antibody, FITC conjugated



PACO61977

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

RNA-binding protein that mediates pre-mRNA alternative splicing regulation. Plays a role in the regulation of cell proliferation, differentiation and migration. Positive regulator of EPO-dependent erythropoiesis. Participates in cell differentiation regulation by repressing tissue-specific exons. Promotes FAS exon 6 skipping. Binds RNA, preferentially to both poly(G) and poly(U).

Gene ID:

PTBP3

Uniprot

O95758

Synonyms:

Polypyrimidine tract-binding protein 3 (Regulator of differentiation 1) (Rod1), PTBP3, ROD1

Immunogen:

Peptide sequence from Human Polypyrimidine tract-binding protein 3 protein (37-55AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A