

# GLU1 Antibody, Biotin conjugated



PACO62006

---

## Product Information

**Size:**

50ul

**Reactivity:**

Arabidopsis thaliana

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in glutamate biosynthesis in leaf. Required for the reassimilation of ammonium ions generated during photorespiration.

**Gene ID:**

GLU1

**Uniprot**

Q9ZLN7

**Synonyms:**

Ferredoxin-dependent glutamate synthase 1, chloroplastic/mitochondrial (EC 1.4.7.1) (Fd-GOGAT 1), GLU1, GLS1

**Immunogen:**

Recombinant Arabidopsis thaliana Ferredoxin-dependent glutamate synthase 1, chloroplastic/mitochondrial protein (106-505AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A