

HMGXB4 Antibody



PACO62103

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Negatively regulates Wnt/beta-catenin signaling during development.

Gene ID:

HMGXB4

Uniprot

Q9UGU5

Synonyms:

HMG domain-containing protein 4 (HMG box-containing protein 4) (High mobility group protein 2-like 1) (Protein HMGBCG), HMGXB4, HMG2L1 HMGBCG


Immunogen:

Peptide sequence from Human HMG domain-containing protein 4 protein (489-509AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: Jurkat whole cell lysate, K562 whole cell lysate. All lanes: HMGXB4 antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 66 kDa. Observed band size: 66 kDa.