

PIGU Antibody, Biotin conjugated



PACO62114

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of the GPI transamidase complex. May be involved in the recognition of either the GPI attachment signal or the lipid portion of GPI.

Gene ID:

PIGU

Uniprot

Q9H490

Synonyms:

Phosphatidylinositol glycan anchor biosynthesis class U protein (Cell division cycle protein 91-like 1) (Protein CDC91-like 1) (GPI transamidase component PIG-U), PIGU, CDC91L1

Immunogen:

Peptide sequence from Human Phosphatidylinositol glycan anchor biosynthesis class U protein (28-47AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A