

AXIN2 Antibody, Biotin conjugated



PACO62122

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Inhibitor of the Wnt signaling pathway. Down-regulates beta-catenin. Probably facilitate the phosphorylation of beta-catenin and APC by GSK3B.

Gene ID:

AXIN2

Uniprot

Q9Y2T1

Synonyms:

Axin-2 (Axin-like protein) (Axil) (Axis inhibition protein 2) (Conductin), AXIN2

Immunogen:

Peptide sequence from Human Axin-2 protein (56-73AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A