

NXN Antibody, FITC conjugated



PACO62153

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Functions as a redox-dependent negative regulator of the Wnt signaling pathway, possibly by preventing ubiquitination of DVL3 by the BCR (KLHL12) complex. May also function as a transcriptional regulator act as a regulator of protein phosphatase 2A (PP2A).

Gene ID:

NXN

Uniprot

Q6DKJ4

Synonyms:

Nucleoredoxin (EC 1.8.1.8), NXN, NRX

Immunogen:

Recombinant Human Nucleoredoxin protein (155-435AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A