

# RPL24 Antibody, HRP conjugated



PACO62160

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

cytoplasm, cytosol, cytosolic large ribosomal subunit, extracellular exosome, membrane, cadherin binding, RNA binding, structural constituent of ribosome, assembly of large subunit precursor of preribosome, nuclear-transcribed mRNA catabolic process, nonsense-mediated decay.

**Gene ID:**

RPL24

**Uniprot**

P83731

**Synonyms:**

60S ribosomal protein L24 (60S ribosomal protein L30) (Large ribosomal subunit protein eL24), RPL24

**Immunogen:**

Recombinant Human 60S ribosomal protein L24 protein (1-157AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A