

RPP38 Antibody, Biotin conjugated



PACO62166

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Component of ribonuclease P, a protein complex that generates mature tRNA molecules by cleaving their 5'-ends. RPP38 may associate transiently with RNase P RNA as a factor involved in the transport of H1 RNA to the putative site of its assembly in the cell, the nucleolus.

Gene ID:

RPP38

Uniprot

P78345

Synonyms:

Ribonuclease P protein subunit p38 (RNaseP protein p38) (EC 3.1.26.5), RPP38

Immunogen:

Recombinant Human Ribonuclease P protein subunit p38 protein (210-283AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A