

PACO62171

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Gene ID:

UBC

Uniprot

P0CG48

Synonyms:

Polyubiquitin-C [Cleaved into: Ubiquitin], UBC


Immunogen:

Peptide sequence from Human Polyubiquitin-C protein (19-36AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: HepG2 whole cell lysate, 293T whole cell lysate, HeLa whole cell lysate, MCF-7 whole cell lysate, K562 whole cell lysate, LO2 whole cell lysate, Jurkat whole cell lysate. All lanes: UBC antibody at 1:2000. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 78 kDa. Observed band size: 78 kDa.