

UBC Antibody, Biotin conjugated



PACO62174

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

UBC

Uniprot

P0CG48

Synonyms:

Polyubiquitin-C [Cleaved into: Ubiquitin], UBC

Immunogen:

Peptide sequence from Human Polyubiquitin-C protein (19-36AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A