

B2m Antibody, Biotin conjugated



PACO62246

Product Information

Size:

50ul

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of the class I major histocompatibility complex (MHC). Involved in the presentation of peptide antigens to the immune system.

Gene ID:

B2m

Uniprot

P01887

Synonyms:

Beta-2-microglobulin, B2m

Immunogen:

Recombinant Mouse Beta-2-microglobulin protein (21-119AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A